



Double washbasin in stainless steelwith surg. lever-operated mixer tap

Available on request

SKU: WASHCH-2 Categories: Sinks

Description

The **stainless-steel double scrub sink** is designed to provide **maximum hygiene and reliability** in operating theatres, clinics and sterilisation areas. Its **all-stainless-steel construction** ensures a robust, corrosion-resistant and easy-to-clean surface, suitable for environments where cleanliness is essential.

The **two-station design** allows simultaneous use by multiple staff members, improving workflow efficiency during surgical preparation. Equipped with **clinical lever mixers**, it enables **hands-free operation**, helping to reduce the risk of cross-contamination.

An integrated **soap dispenser** ensures convenient access during handwashing, supporting consistent hygiene practices.

The **deep bowl with rounded edges** and **smooth welded surfaces** prevent residue build-up and make cleaning quick and efficient.

Built to deliver **durability**, **safety and ease of use**, this stainless-steel double scrub sink provides a **professional solution for hospitals**, **surgical units and clinics** where hygiene and performance are critical requirements.