



DOVA – Dental Office Voice Assistant for Windows 10-11

D.O.V.A. is an advanced voice assistant built specifically for dental and medical practices, allowing staff to manage appointments, clinical records and notes completely hands-free.

Estimated shipping in 10 days

SKU: FK811MTR

Categories: All-in-One Medical PCs

Description

D.O.V.A. – Dental Office Voice Assistant

D.O.V.A. – Dental Office Voice Assistant is a professional voice assistant designed to significantly reduce the time spent on clinical documentation through **high-precision voice dictation**.

It is an ideal solution for dental practices and medical environments that manage a high volume of clinical and operational information on a daily basis.

During clinical and organisational activities, D.O.V.A. allows spoken information to be converted into digital text **in real time**. This improves workflow continuity and the overall organisation of clinical documentation.

Thanks to its **plug-and-play design**, D.O.V.A. interfaces with applications running on the PC. It enables voice dictation to be used directly within the software already in use, without altering existing workflows or requiring complex configuration.

Key features and benefits

High-precision voice dictation

Designed for professional healthcare environments, D.O.V.A. provides accurate and reliable voice dictation. It is suitable for clinical records, operative notes and medical documentation.

Plug-and-play operation

D.O.V.A. is ready to use and works across applications installed on the PC. It allows voice dictation within existing software without additional installations or workflow changes.



Wireless microphones for flexible use

The wireless microphones allow operation across different rooms and workstations. They adapt easily to the daily routines of clinicians and support staff.

Reduced documentation time

Immediate digitisation of spoken information helps reduce the time dedicated to documentation. Information remains available and up to date without delays.

Support for dental assistants and clinical teams

D.O.V.A. supports the work of dental assistants and clinical staff. It improves internal organisation and facilitates collaboration within the practice.

Workflow continuity

The voice assistant supports daily activities while minimising interruptions. This helps professionals remain focused on clinical procedures.

Versatile application

Suitable for a wide range of healthcare settings, D.O.V.A. fits naturally into existing workflows without the need for structural changes.

Technical specifications

- **Technology:** Professional voice assistant
- **Interaction:** High-precision voice dictation
- **Microphones:** Two wireless microphones included
- **Foot pedals:** Three configurable pedals
- **Connection:** Wireless
- **Compatibility:** PC
- **Intended use:** Dental practices and medical environments