



Dry air compressor, 1 unit, 70 l/min, without dryer

Tecnomed Silent Line dry air compressor for 1 dental unit, 25l tank, air output 70 l/min, without dryer

Available

SKU: 100TOV

Categories: Without air dryer

Description

Medical Dry-air compressor 1 dental unit 70l/min.

Tecnomed Silent Compressors – Unmatched Silence and Performance for Your Practice

The Tecnomed Silent Dry Compressor is the perfect choice for dental practices and clinics that require a silent solution for 1 dental units. Engineered to deliver a serene working environment, high efficiency, and long-lasting durability, this compressor is designed with your practice's needs in mind.

Key Features:

- **Exceptional Quietness:** The Tecnomed Silent Compressor operates with minimal noise, creating a more comfortable work environment.
- **Optimal Performance:** Equipped with a durable copper motor and oil-free piston, this compressor ensures high efficiency and reduced heat dissipation.
- **Superior Durability:** With a hard anodized process and anti-corrosion internal coating, the compressor is built to last.

Why Choose Tecnomed Silent Compressors?

Tecnomed Silent Compressors are ideal for those seeking top-notch quietness and performance. Featuring advanced technology, these compressors ensure a peaceful and comfortable working environment while maintaining high standards of efficiency and durability. They provide purified air and long operational life, making them the optimal choice for dental practices and clinics that demand quality and reliability without compromise.

Technical features Medical Dry-air compressor 1 dental unit 70l/min.



Dental units served	1
Air flow70	liters/minutes
Aspirated volume	115 liters
Tank capability	25 liters
Oil less	•
Air dryer	No
Own cabinet	
Twin cylinder stages	1
Power	550 Watt
Noise level	52 db
Power supply	220/240V-50/60 Hz
Weight	24 kg
Dimensions (W x L x H)	400 x 400 x 530 mm
Warranty	24 months
Certification	CE