



## Dry air compressor, 2 units, 140 l/min without dryer

Tecnomed Silent Line dry air compressor for 2 treatment units, 38l tank, air output 140 l/min, without dryer

**Estimated shipping in 60 days**

**SKU:** 240TOV

**Categories:** Without Dryer

### Description

---

#### Medical Air Compressor 2 Units – Silent, Efficient and Reliable for Dental Clinics

The **Tecnomed Silent** medical air compressor 2 units is designed to supply clean, dry and safe compressed air for two dental units. It ensures consistent performance and a quiet working environment, perfectly suited to the needs of modern dental practices and medical facilities.

#### Silent and Comfortable Operation

The soundproofed system guarantees exceptionally quiet performance, reducing noise levels and improving the acoustic comfort of the workspace. This creates a calm and professional environment, enhancing operator focus and ensuring a more relaxed experience for patients.

#### Clean Air and Oil-Free Technology

The **completely oil-free head**, compliant with medical and dental standards, delivers pure and uncontaminated compressed air. This feature reduces maintenance and ensures maximum hygiene and reliability during all clinical applications.

#### Durability and Long-Term Resistance

The tank undergoes **hard anodising** and features an **internal anti-corrosion coating**, which increase strength and service life, even under intensive operating conditions.

#### Reliable Choice for Modern Practices



Silent, efficient and long-lasting, the **Tecnomed Silent** system provides a steady flow of clean air for dental surgeries and clinics that demand consistent performance and full compliance with medical quality standards.

**Technical features Medical Dry-air compressor 2 dental units.**

Dental units served	2
Air flow	140 liters/minutes
Aspirated volume	230 liters
Tank capability	38 liters
Oil less	•
Air dryer	No
Own cabinet	
Twin cylinder stages	1
Power	1100 Watt
Noise level	52 db
Operating Pressure	5,5-8 Bar
Power supply	220/240V-50/60 Hz
Weight	34 kg
Dimensions (W x L x H)	450 x 450 x 750 mm
Warranty	24 months
Certification	CE