



Dry air compressor, 3 units, 195 l/min, with dryer

Tecnomed Silent Line dry air compressor for 3 treatment units, 60l tank, air output 195 l/min, with dryer

Available on request

SKU: 360TOVEN

Categories: With Dryer

Description

Medical Air Compressor 3 Units with Dryer – Quiet, Efficient and Reliable for Dental Clinics

The Tecnomed Silent medical air compressor 3 units with dryer delivers clean, dry and reliable compressed air for 3 dental units. Designed for professional environments, it combines silent operation, steady performance and long-term dependability — making it ideal for modern dental surgeries and medical practices.

Whisper-Quiet Operation

The advanced soundproofing system ensures extremely quiet performance, creating a calm and comfortable working atmosphere. Lower noise levels help maintain staff focus and provide patients with a more relaxing experience during treatment.

Oil-Free and Hygienic Air Supply

Equipped with a fully oil-free head, the compressor provides pure, uncontaminated air that meets strict medical and dental standards. This technology eliminates the risk of oil residues, minimises maintenance and guarantees optimal hygiene for all clinical applications.

Integrated Dryer System

The integrated dryer features a high-efficiency multi-stage filtration process that removes humidity and residual particles from the compressed air. This ensures a constant flow of dry, clean air that protects turbines, handpieces and internal components from wear, corrosion and moisture-related damage — maintaining consistent performance over time.



Built to Last

A hard-anodised tank with an internal anti-corrosion coating offers maximum protection against wear and moisture. These features ensure consistent performance even under intensive daily use.

Dependable Performance for Modern Practices

Quiet, efficient and easy to maintain, the Tecnomed Silent compressor with dryer ensures a constant supply of high-quality air — the reliable choice for dental professionals who demand safety, precision and consistency in every procedure.

Technical features Medical Dry air compressor 3 dental units with dryer.

Dental units served	3
Air flow	195 liters/minutes
Aspirated volume	304 liters
Tank capability	60 liters
Oil less	•
Air dryer	•
Own cabinet	
Twin cylinder stages	2
Power	1500 Watt
Noise level	53 db
Operating Pressure	5,5-8 Bar
Power supply	220/240V-50/60 Hz
Weight	37 kg
Dimensions (W x L x H)	715 x 500 x 735 mm
Warranty	24 months
Certification	CE