



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/dry-air-compressor-3-units-205-l-min-without-dryer/>



## Dry air compressor, 3 units, 205 l/min without dryer

Tecnomed Silent Line dry air compressor for 3 treatment units, 60l tank, air output 205 l/min, without dryer

**1.886,51 €**

**Available**

**SKU:** 360TOV

**Categories:** Without air dryer

### Description

Medical Dry-air compressor 3 dental units.

### Tecnomed Silent Compressors – Unmatched Silence and Performance for Your Practice

The Tecnomed Silent Dry Compressor is the perfect choice for dental practices and clinics that require a silent solution for 3 dental units. Engineered to deliver a serene working environment, high efficiency, and long-lasting durability, this compressor is designed with your practice's needs in mind.

#### Key Features:

- **Exceptional Quietness:** The Tecnomed Silent Compressor operates with minimal noise, creating a more comfortable work environment.
- **Optimal Performance:** Equipped with a durable copper motor and oil-free piston, this compressor ensures high efficiency and reduced heat dissipation.
- **Superior Durability:** With a hard anodized process and anti-corrosion internal coating, the compressor is built to last.

### Why Choose Tecnomed Silent Compressors?

Tecnomed Silent Compressors are ideal for those seeking top-notch quietness and performance. Featuring advanced technology, these compressors ensure a peaceful and comfortable working environment while maintaining high standards of efficiency and durability. They provide purified air and long operational life, making them the optimal choice for dental practices and clinics that demand quality and reliability without compromise.

#### Technical features Medical Dry-air compressor 3 dental units.

Dental units served	3
Air flow	195 liters/minutes
Aspirated volume	304 liters

Tank capability	60 liters
Oil less	•
Air dryer	No
Own cabinet	
Twin cylinder stages	2
Power	1500 Watt
Noise level	52 db
Power supply	220/240V-50/60 Hz
Weight	51 kg
Dimensions (W x L x H)	715 x 367 x 735 mm
Warranty	24 months
Certification	CE