



# Dry air compressor for 1+1 dental units, 95l/min, with dryer

Tecnomed Silent Line dry air compressor for 1+1 dental units, 38l tank, air output 95 l/min, with dryer

#### **Available**

SKU: 150TOVEN

Categories: With Dryer

# Description

# Medical Air Compressor 1+1 Units with Dryer – Quiet, Efficient and Reliable for Dental Clinics

The Tecnomed Silent medical air compressor 1+1 units with dryer delivers clean, dry and reliable compressed air for one or two dental units. Designed for professional environments, it combines silent operation, steady performance and long-term dependability — making it ideal for modern dental surgeries and medical practices.

#### **Whisper-Quiet Operation**

The advanced soundproofing system ensures extremely quiet performance, creating a calm and comfortable working atmosphere. Lower noise levels help maintain staff focus and provide patients with a more relaxing experience during treatment.

#### Oil-Free and Hygienic Air Supply

Equipped with a fully oil-free head, the compressor provides pure, uncontaminated air that meets strict medical and dental standards. This technology eliminates the risk of oil residues, minimises maintenance and guarantees optimal hygiene for all clinical applications.

# **Integrated Dryer System**

The integrated dryer features a high-efficiency multi-stage filtration process that removes humidity and residual particles from the compressed air. This ensures a constant flow of dry, clean air that protects turbines, handpieces and internal components from wear, corrosion and moisture-related damage —



maintaining consistent performance over time.

#### **Built to Last**

A hard-anodised tank with an internal anti-corrosion coating offers maximum protection against wear and moisture. These features ensure consistent performance even under intensive daily use.

## **Dependable Performance for Modern Practices**

Quiet, efficient and easy to maintain, the Tecnomed Silent compressor with dryer ensures a constant supply of high-quality air — the reliable choice for dental professionals who demand safety, precision and consistency in every procedure.

## Technical features Medical dry-air compressor 1+1 dental units with dryer

Dental units served	1+1
Air flow	95 liters/minutes
Aspirated volume	152 liters
Tank capability	38 liters
Oil less	•
Air dryer	•
Own cabinet	
Twin cylinder stages	1
Power	750 Watt
Noise level	53 db
Power supply	220/240V-50/60 Hz
Weight	37 kg
Dimensions (W x L x H)	430 x 400 x 660 mm
Warranty	24 months
Certification	CE