

## Dry air compressor for 5 dental unit Tecnomed Silent Line, tank 100 LT



Tecnomed Silent Line dry air compressor for 5 treatment units, 100l tank, air output 312 l/min, with dryer

**Shipping expected in approximately 60 days**

**SKU:** 720TAIN

**Categories:** With Dryer

### Description

---

#### Medical Air Compressor 5Units with Dryer – Quiet, Efficient and Reliable for Dental Clinics

The Tecnomed Silent medical air compressor 5 units with dryer delivers clean, dry and reliable compressed air for up to 5 dental units. Designed for professional environments, it combines silent operation, steady performance and long-term dependability — making it ideal for modern dental surgeries and medical practices.

#### Whisper-Quiet Operation

The advanced soundproofing system ensures extremely quiet performance, creating a calm and comfortable working atmosphere. Lower noise levels help maintain staff focus and provide patients with a more relaxing experience during treatment.

#### Oil-Free and Hygienic Air Supply

Equipped with a fully oil-free head, the compressor provides pure, uncontaminated air that meets strict medical and dental standards. This technology eliminates the risk of oil residues, minimises maintenance and guarantees optimal hygiene for all clinical applications.

#### Integrated Dryer System

The integrated dryer features a high-efficiency multi-stage filtration process that removes humidity and residual particles from the compressed air. This ensures a constant flow of dry, clean air that protects turbines, handpieces and internal components from wear, corrosion and moisture-related damage — maintaining consistent performance over time.



### **Built to Last**

A hard-anodised tank with an internal anti-corrosion coating offers maximum protection against wear and moisture. These features ensure consistent performance even under intensive daily use.

### **Dependable Performance for Modern Practices**

Quiet, efficient and easy to maintain, the Tecnomed Silent compressor with dryer ensures a constant supply of high-quality air — the reliable choice for dental professionals who demand safety, precision and consistency in every procedure.

### **Technical features Medical Dry-air compressor 5 dental unit with dryer.**

Dental units served	5
Air flow	312 liters/minutes
Aspirated volume	720 liters
Tank capability	100 liters
Air Dryer	Yes
Power	3000 Watt
Noise level	71 db
Operating Pressure	5,5-8 Bar
Weight	94 kg