



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/easy-tray-nano-cart-for-instruments-with-3-trays-465x315-mm/>



## Easy Tray Nano cart for instruments, with 3 trays 465x315 mm

**1.061,73 €**

**Available on request**

**SKU:** OF1003BIN.01

**Categories:** Multifunctional cart

### Description

---

Easy Tray Nano 3 trays cart for equipment is the best solution for host you medical instruments in a orderly manner.

**Versatile and modular:** Made entirely of electro galvanized white enamelled sheet, this cart is modular and configurable with different accessories according to your needs.

**No cable in your way:** you'll get a cleaner, more sterile environment thanks to hidden electrical connections and a concealed power supply.

**One cart for all your equipment:** 3 sterilizable and removable trays for host complementary surgical equipment such as scalers, electrosurgery units, piezoelectric surgical devices, implant micromotors, endodontic micromotors, lasers, surgical Stryker equipment, whitening devices, multi-parameter monitors. A single, integrated, height-adjustable infusion set stand can hold up to three bags of saline solution to ensure easy operation.

Technical feature **Easy Tray Nano 3 trays cart:**

- height 85,5 cm;
- 3 stainless steel tray (46,5x31,5x2 cm), removable and sterilizable;
- Centralized power system for equipment.