



## Easy Tray Nano surgical equipment trolley white, 3 removable tray

**Available**

**SKU:** OF1003BIN.01

**Categories:** Equipment Carts

### Description

---

Easy Tray Nano White is a compact powered medical trolley manufactured in white powder-coated steel, designed for the organised placement and power supply of electromedical, dental and surgical equipment in clinical and operating environments.

Developed for dental practices, veterinary clinics, outpatient facilities and operating theatres, it allows electrical devices to be used directly on the trolley, ensuring order, electrical safety and hygiene even in limited spaces. The white finish integrates naturally into modern healthcare settings while delivering the same functional performance as the stainless steel version.

Thanks to the integrated power supply, Easy Tray Nano White reduces cable clutter and enhances operational safety, simplifying the integration of equipment directly on the trolley. Concealed electrical cable management helps keep the work area clean, orderly and easy to disinfect, improving workflow efficiency and ergonomics during clinical procedures.

The compact and stable structure makes Easy Tray Nano White suitable for intensive clinical use, providing long-term reliability and ease of use without unnecessary complexity.

---

### Key features

- Powered medical trolley for clinical and electromedical equipment
- White powder-coated steel structure, robust and easy to clean and disinfect
- 3 removable and sterilizable trays (465 x 315 x 20 mm)



- Concealed electrical cable management for improved safety and organisation
- Compact Nano design, ideal for practices with optimised working spaces
- Medical-grade castors with brakes for stability during procedures
- Suitable as a medical trolley for clinics, dental practices and surgical environments

Overall dimensions:

Height 87.5cm x Width 60cm x Depth 50cm

---

### **Optional accessories (available separately)**

- Additional tray 415 x 300 x 20 mm (ref. DE1003X.02)
- Extension for Easy Tray Nano White with 2 removable and sterilizable trays 465 x 315 x 20 mm (ref. OF1003BIN.02)