



## Easy Tray Nano surgical equipment trolley white, 3 removable tray

**Available on request**

**SKU:** OF1003BIN.01

**Categories:** Equipment Carts

### Description

---

Easy Tray Nano surgical equipment trolley with 3 trays is the best solution for host you medical instruments in a orderly manner.

**Versatile and modular:** This cart is modular and configurable with different accessories according to your needs.

**No cable in your way:** you'll get a cleaner, more sterile environment thanks to hidden electrical connections and a concealed power supply.

**One cart for all your equipment:** 3 sterilizable and removable trays for host complementary surgical equipment such as scalers, electrosurgery units, piezoelectric surgical devices, implant micromotors, endodontic micromotors, lasers, surgical Stryker equipment, whitening devices, multi-parameter monitors. A single, integrated, height-adjustable infusion set stand can hold up to three bags of saline solution to ensure easy operation.

**Fields of use:** ideal for medical, surgical and dental settings, providing a safe and well-organised way to manage clinical equipment.

#### **Technical feature Easy Tray Nano surgical equipment trolley:**

- height 85,5 cm;
- 3 stainless steel tray (46,5x31,5x2 cm), removable and sterilizable;
- Centralized power system for equipment.