



## Easy Tray Wall Electrified Medical Equipment Station

Electrified wall-mounted station for medical equipment with articulated arm and 3 removable stainless steel trays, designed to support and power devices in dental practices and clinical environments.

**Available on request**

**SKU:** OF1003BP

**Categories:** Equipment Carts

### Description

---

#### Easy Tray Wall electrified medical equipment wall station

---

**Easy Tray Wall** is an **electrified wall-mounted station for medical equipment**, designed to support and power devices in dental practices, outpatient settings and clinical environments.

The integrated electrical system allows direct connection of medical devices, ensuring a clean, organized and efficient workspace during clinical procedures.

The **double-jointed articulated arm** enables precise and flexible positioning of equipment, improving ergonomics and accessibility during daily operations.

The system includes **3 removable trays in medical-grade stainless steel**, designed to organize instruments and devices efficiently while ensuring durability, hygiene and easy surface disinfection.

The **wall-mounted configuration** optimizes available space, keeping equipment accessible without occupying floor area and improving overall workflow organization.

This solution follows the design logic of equipment support systems used in operating rooms, adapted to the needs of dental practices and clinical environments.

**Easy Tray Wall** provides a compact and efficient solution for medical equipment management, enhancing ergonomics, workflow efficiency and clinical organization.

### Key features

- **wall-mounted support for medical equipment**
- **integrated electrified structure with power outlets**



- **double-jointed articulated arm**
- **3 removable medical-grade stainless steel trays**
- surfaces designed for **easy cleaning and disinfection**
- durable coated structure for clinical environments
- ideal for **dental practices, clinics and outpatient settings**