



## Easy Tray white surgical equipment trolley, without trays

Surgical equipment trolley Easy Tray white, without trays / shelves

**Available**

**SKU:** DE1003BI.00

**Categories:** Equipment Carts

### Description

---

Easy Tray White is a powered medical trolley made of white powder-coated steel, designed for powering and organising electromedical, dental and surgical equipment in high-intensity clinical and operating environments. The white finish blends naturally into modern healthcare settings while delivering the same functional performance as the stainless steel version.

Easy Tray White is suitable for dental practices, outpatient clinics, operating theatres and hospital environments where an electrified mobile base is required to support multiple devices during clinical procedures. Thanks to the integrated power supply, equipment can be powered directly from the trolley, reducing cable clutter and enhancing operational safety.

Concealed cable management keeps the workspace clean, organised and easy to sanitise. The solid and stable construction makes Easy Tray White suitable for intensive clinical use, ensuring long-term reliability without unnecessary complexity.

Easy Tray White acts as a **modular platform**, configurable with a wide range of optional trays and accessories to adapt the trolley precisely to professional workflow requirements.

***This version of Easy Tray White is supplied without trays, allowing users to freely configure the trolley by selecting the trays that best suit their specific operational needs. The structure is preconfigured to accommodate Easy Tray pull-out trays, enabling a fully modular and customised workstation setup.***

### Key features

- Powered medical trolley for clinical and electromedical equipment



- White powder-coated steel structure, durable and easy to sanitise
- Designed to accommodate Easy Tray pull-out trays (not included)
- Concealed electrical connections for improved safety and order
- Robust and stable design for intensive clinical use
- Medical-grade castors with brakes for procedural stability

### **Dimensions**

- Height: 145cm (215cm with IV pole fully extended)
- Width: 60cm
- Depth: 50cm

### **Optional accessories (to be ordered separately)**

- Easy Tray pull-out tray 415×300×20mm (**code DE1003X.02**)
- Easy Tray pull-out tray 465×315×20mm (**code DE1003X.01**)
- 5-tray instrument kit 415×300×20mm (**code OF1003X.2**)
- 5-tray instrument kit 465×315×20mm (**code OF1003X.1**)
- Monitor support with VESA 75/100 mount for medical monitors (**ref. SVET9010L**)
- Handpiece arm for 5-tray Easy Tray (**ref. DE1003X.5M**)
- Implantology handpiece holder (**ref. DE1003X.1**)
- Piezosurgery handpiece holder (**ref. DE1003X.2**)



- Electrosurgery handpiece holder (ref. **DE1003X.3**)
- OKKIO medical camera arm (ref. **DE1003X.6K**)