



Easy Tray inox surgical equipment trolley, without trays

Surgical powered trolley structure made of stainless steel, without trays, equipped with handle, IV pole and 4 swivel castors with brake.

Available

SKU: DE1003X.00

Categories: Equipment Carts

Description

Easy Tray Inox is a powered surgical equipment trolley (version without trays) made of stainless steel, designed for organising and powering electromedical, dental and surgical equipment in high-intensity clinical and operating environments. The stainless steel construction ensures durability, hygiene and easy sanitisation, making the trolley suitable for intensive professional use.

Easy Tray Inox is suitable for dental practices, outpatient clinics, operating theatres and hospital environments where an electrified mobile base is required to support multiple devices during clinical procedures. Thanks to the integrated power supply, equipment can be powered directly from the trolley, reducing cable clutter and improving operational safety.

Concealed cable management helps keep the working area clean, organised and easy to sanitise. The solid and stable construction ensures long-term reliability without adding unnecessary complexity to daily use.

Easy Tray Inox functions as a **modular platform**, configurable with optional trays and accessories, allowing the trolley to be adapted precisely to professional workflow needs.

This version of Easy Tray Inox is supplied without trays, allowing professionals to freely configure the workstation by selecting the trays that best match their specific operational requirements. The trolley structure is designed to accommodate Easy Tray pull-out trays, enabling a modular and fully customisable setup.

Key features

- Powered medical trolley for clinical and electromedical equipment



- Stainless steel structure, durable and easy to sanitise
- Designed to accommodate Easy Tray pull-out trays (not included)
- Concealed electrical connections for improved safety and order
- Robust and stable design for intensive clinical use
- Medical-grade castors with brakes for stability during procedures
- Suitable as a medical trolley for dental, clinical and surgical environments

Dimensions

- Height: 145cm (215cm with IV pole fully extended)
- Width: 60cm
- Depth: 50cm

Optional accessories (to be ordered separately)

- Easy Tray pull-out tray 415×300×20mm (**code DE1003X.02**)
- Easy Tray pull-out tray 465×315×20mm (**code DE1003X.01**)
- 5-tray instrument kit 415×300×20mm (**code OF1003X.2**)
- 5-tray instrument kit 465×315×20mm (**code OF1003X.1**)
- Monitor support with VESA 75/100 mount for medical monitors (**ref. SVET9010L**)
- Handpiece arm for 5-tray Easy Tray (**ref. DE1003X.5M**)
- Implantology handpiece holder (**ref. DE1003X.1**)



- Piezosurgery handpiece holder (**ref. DE1003X.2**)
- Electrosurgery handpiece holder (**ref. DE1003X.3**)