



Printed on 01/06/2025

<https://shop.tecnomeditalia.com/en/product/endo-tip-for-scaler-compatible-satelec-120-dte-e1d-nsk-e11/>



Endo tip for scaler, compatible Satelec 120, DTE E1D, NSK E11

Endo tip 120° for root canal cleaning of molars. It accommodates all flexible files. Compatible Satelec 120, DTE E1D, NSK E11

29,14 €

Available

SKU: PG115

Categories: Endodontic inserts

Description

Endo tip for scaler compatible Satelec-120 DTE-E1D NSK-E11.

Tecnomed Italia endodontic inserts are the perfect complement for rotary instruments in the preparation of endodontic cavities. Extremely thin in shape, they perfectly respect canalar anatomy and provide easy access to pulp chambers.

They are also ideal for removing filling residues and preparing and irrigating the canalar area. The files are extremely flexible and strong, and adapt to all canalar shapes, reducing the risk of damage.

This makes them ideal for shaping, cleaning and disinfecting endodontic cavities.