



Endo tip for scaler, compatible Satelec 120, DTE E1D, NSK E11

Endo tip 120° for root canal cleaning of molars. It accommodates all flexible files. Compatible Satelec 120, DTE E1D, NSK E11

Available

SKU: PG115

Categories: Satelec - DTE - NSK

Description

Endo tip for scaler compatible with Satelec 120, DTE E1D and NSK E11, for root canal cleaning of molars.

Specifications

- Compatible with Satelec 120, DTE E1D and NSK E11
- Suitable for use with compatible Satelec, DTE and NSK piezoelectric scalars
- For professional dental use