



Printed on 01/06/2025

<https://shop.tecnomeditalia.com/en/product/endo-tip-for-scaler-compatible-satelec-90-dte-e2d-nsk-e12/>



## Endo tip for scaler, compatible Satelec 90, DTE E2D, NSK E12

Endo tip for root canal cleaning of molars. It accommodates all flexible files. Compatible Satelec 90, DTE E2D, NSK E12

**22,72 €**

**Available**

**SKU:** PG116

**Categories:** Endodontic inserts

### Description

---

Endo tip for scaler compatible Satelec-90 DTE-E2D NSK-E12.

Tecnomed Italia endodontic inserts are the perfect complement for rotary instruments in the preparation of endodontic cavities. Extremely thin in shape, they perfectly respect canalar anatomy and provide easy access to pulp chambers.

They are also ideal for removing filling residues and preparing and irrigating the canalar area.

The files are extremely flexible and strong, and adapt to all canalar shapes, reducing the risk of damage. This makes them ideal for shaping, cleaning and disinfecting endodontic cavities.