


Endo tip for scaler, EMS compatible 120, Woodpecker E1



Endo tip 120° for root canal cleaning of molars. It accommodates all flexible files. EMS compatible 120, Woodpecker E1

Available

SKU: PG066

Categories: EMS - Woodpecker

Description

Endo tip for scaler compatible EMS 120 Woodpecker E1.

Tecnomed Italia endodontic inserts are the perfect complement for rotary instruments in the preparation of endodontic cavities. Extremely thin in shape, they perfectly respect canal anatomy and provide easy access to pulp chambers.

They are also ideal for removing filling residues and preparing and irrigating the canal area. The files are extremely flexible and strong, and adapt to all canal shapes, reducing the risk of damage.

This makes them ideal for shaping, cleaning and disinfecting endodontic cavities.

Compatibility Endo tip for scaler: EMS 120, Woodpecker E1.