


Endo tip for scaler, EMS compatible 120, Woodpecker E1



Endo tip 120° for root canal cleaning of molars. It accommodates all flexible files. EMS compatible 120, Woodpecker E1

Available

SKU: PG066

Categories: EMS - Woodpecker

Description

The 120° endo scaler tip, EMS and Woodpecker E1 compatible is designed for root canal cleaning in molars, providing an angled shape suitable for access to posterior areas.

It is designed to hold flexible files, supporting endodontic cleaning procedures inside the root canal.

Specifications

- Compatible with EMS 120 and Woodpecker E1
- Suitable for use with compatible piezoelectric scalars
- For professional dental use