



Printed on 01/06/2025

<https://shop.tecnomeditalia.com/en/product/endo-tip-for-scaler-ems-compatible-90-woodpecker-e2/>



Endo tip for scaler, EMS compatible 90, Woodpecker E2

Endo tip for root canal cleaning of molars. It accommodates all flexible files. EMS 90, Woodpecker E2 compatible

22,72 €

Available

SKU: PG069

Categories: Endodontic inserts

Description

Endo tip for scaler compatible EMS 90 Woodpecker E2.

Tecnomed Italia endodontic inserts are the perfect complement for rotary instruments in the preparation of endodontic cavities.

Extremely thin in shape, they perfectly respect canal anatomy and provide easy access to pulp chambers. They are also ideal for removing filling residues and preparing and irrigating the canal area.

The files are extremely flexible and strong, and adapt to all canal shapes, reducing the risk of damage. This makes them ideal for shaping, cleaning and disinfecting endodontic cavities.