



Endo tip for scaler, EMS compatible 90, Woodpecker E2

Endo tip for root canal cleaning of molars. It accommodates all flexible files. EMS 90, Woodpecker E2 compatible

Shipping expected in approximately 40 days

SKU: PG069

Categories: EMS - Woodpecker

Description

The endo tip for scaler compatible with EMS 90 and Woodpecker E2 is designed for root canal cleaning in molars.

The tip allows the use of flexible files, supporting professional endodontic procedures for cleaning and preparing the root canal space.

Specifications

- Compatible with EMS 90 and Woodpecker E2
- Suitable for use with EMS and Woodpecker compatible piezoelectric scalars
- For professional dental use