



Endo tip for scaler, EMS compatibleH, Woodpecker E3

Endo tip for lateral thermal condensation of gutta-percha. EMS H, Woodpecker E3 compatible

Available

SKU: PG054

Categories: Endodontic inserts

Description

Endo tip for scaler compatible EMS H Woodpecker E3.

Tecnomed Italia endodontic inserts are the perfect complement for rotary instruments in the preparation of endodontic cavities. Extremely thin in shape, they perfectly respect canalar anatomy and provide easy access to pulp chambers.

They are also ideal for removing filling residues and preparing and irrigating the canalar area. The files are extremely flexible and strong, and adapt to all canalar shapes, reducing the risk of damage.

This makes them ideal for shaping, cleaning and disinfecting endodontic cavities.