



Printed on 01/06/2025

<https://shop.tecnomeditalia.com/en/product/endo-tip-for-scaler-ems-compatibleh-woodpecker-e3/>



## Endo tip for scaler, EMS compatibleH, Woodpecker E3

Endo tip for lateral thermal condensation of gutta-percha. EMS H, Woodpecker E3 compatible

**24,11 €**

**Available**

**SKU:** PG054

**Categories:** Endodontic inserts

### Description

---

Endo tip for scaler compatible EMS H Woodpecker E3.

Tecnomed Italia endodontic inserts are the perfect complement for rotary instruments in the preparation of endodontic cavities. Extremely thin in shape, they perfectly respect canal anatomy and provide easy access to pulp chambers.

They are also ideal for removing filling residues and preparing and irrigating the canal area. The files are extremely flexible and strong, and adapt to all canal shapes, reducing the risk of damage.

This makes them ideal for shaping, cleaning and disinfecting endodontic cavities.