



Endo tip for scaler, Mectron D3 compatible

Endo tip for lateral thermal condensation of gutta-percha. Mectron D3 compatible

Available

SKU: PG071M

Categories: Endodontic inserts

Description

Endo tip for scaler Mectron D3 compatible.

Tecnomed Italia endodontic inserts are the perfect complement for rotary instruments in the preparation of endodontic cavities.

Extremely thin in shape, they perfectly respect canal anatomy and provide easy access to pulp chambers. They are also ideal for removing filling residues and preparing and irrigating the canal area.

The files are extremely flexible and strong, and adapt to all canal shapes, reducing the risk of damage. This makes them ideal for shaping, cleaning and disinfecting endodontic cavities.