



ERGO A seat, 5003 midnight blue, aluminium base

ERGO A seat, colour 5003 midnight blue, aluminium base, height adjustment and backrest rotation

Available from 11/05/2026

SKU: DE8.06T5003

Categories: Over Line

Description

ERGO A seat 5003 midnight blue.

An appropriate stool is an important investment for the comfort and health of the dental professional. Tecnomed has a vast experience in the design and production of ergonomic dental stools. Tecnomed stools help you maintain a correct posture with natural curvature of the spine, uniform distribution of body weight and no painful compressions.

ERGO A stools are designed to offer superior comfort for dentists and patients at all times.

ERGO A stools are upholstered in top quality imitation leather and are free from stitching for maximum hygiene and easy cleaning.

The upholstery material is fireproof according to ISO 6875 and UNI 4818.74. This top quality upholstery helps ensure a soft, anatomical seating position with uniform weight distribution and no painful compressions.

The “U” shaped, rotating backrest on the ERGO A models provides good lumbar support and a practical rest for the user’s arm. The backrest can be reversed, right to left, making these stools the perfect choice even for left handed dentists.

ERGO A seat 5003 midnight blue.