



Printed on 07/06/2025

<https://shop.tecnomeditalia.com/en/product/ergo-a-seat-5012-periwinkle-blue-aluminium-base/>



ERGO A seat, 5012 periwinkle blue, aluminium base

ERGO A seat, colour 5012 periwinkle blue, aluminium base, height adjustment and backrest rotation

870,10 €

Available

SKU: DE8.06T3491

Categories: Comfort

Description

ERGO A seat 5012 periwinkle blue.

An appropriate stool is an important investment for the comfort and health of the dental professional. Tecnomed has a vast experience in the design and production of ergonomic dental stools. Tecnomed stools help you maintain a correct posture with natural curvature of the spine, uniform distribution of body weight and no painful compressions.

ERGO A stools are designed to offer superior comfort for dentists and patients at all times.

ERGO A stools are upholstered in top quality imitation leather and are free from stitching for maximum hygiene and easy cleaning.

The upholstery material is fireproof according to ISO 6875 and UNI 4818.74. This top quality upholstery helps ensure a soft, anatomical seating position with uniform weight distribution and no painful compressions.

The "U" shaped, rotating backrest on the ERGO A models provides good lumbar support and a practical rest for the user's arm. The backrest can be reversed, right to left, making these stools the perfect choice even for left handed dentists.

ERGO A seat 5012 periwinkle blue.