



ERGO A seat, 6019 apple green, aluminium base

ERGO A seat, colour 6019 apple green, aluminium base, height adjustment and backrest rotation

Available on request

SKU: DE8.06T5515

Categories: Comfort

Description

ERGO A seat 6019 apple green.

An appropriate stool is an important investment for the comfort and health of the dental professional. Tecnomed has a vast experience in the design and production of ergonomic dental stools. Tecnomed stools help you maintain a correct posture with natural curvature of the spine, uniform distribution of body weight and no painful compressions.

ERGO A stools are designed to offer superior comfort for dentists and patients at all times.

ERGO A stools are upholstered in top quality imitation leather and are free from stitching for maximum hygiene and easy cleaning.

The upholstery material is fireproof according to ISO 6875 and UNI 4818.74. This top quality upholstery helps ensure a soft, anatomical seating position with uniform weight distribution and no painful compressions.

The “U” shaped, rotating backrest on the ERGO A models provides good lumbar support and a practical rest for the user’s arm. The backrest can be reversed, right to left, making these stools the perfect choice even for left handed dentists.

ERGO A seat 6019 apple green.