



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/falkoled-light-5700k-24000-lux-jack-connection-for-okkio-hdmi/>



FALKOLED Light 5700K°, 24000 LUX, jack connection, for OKKIO HDMI

LED illuminator designed for OKKIO HDMI cameras, round spotlight of 24,000 lux, 5700 K° color temperature, brightness adjustment via plug-in touchpad, powered by 12 VDC.

472,91 €

Available

SKU: IL901H

Categories: Medical Cameras

Description

FALKOLED Illuminator for OKKIO HDMI.

Cameras Lighting quality is essential for capturing high-quality photos and videos. The IL901H, specifically designed for OKKIO HDMI cameras, offers a round spotlight of 24,000 lux and a color temperature of 5700 K°, providing clear, shadow-free visibility, ideal for recording in clinical environments. The plug-in touchpad allows easy brightness adjustment, ensuring the light is always suited to the recording conditions.

High-quality lighting is crucial to avoid glare, shadows, and color distortions that could affect image quality. The IL901H ensures stable illumination, which is essential for medical procedures requiring visual precision.

Powered by 12 VDC, this illuminator provides the power necessary for professional, extended use, delivering excellent visual results in any situation.

Technical features:

- Specifically designed for OKKIO HDMI cameras.
- round spotlight with 24,000 lux brightness.
- 5700 K° color temperature for natural, stable light.
- On/off and brightness adjustment via plug-in touchpad.
- Powered by 12 VDC.
- Stable, even light for high-quality videos and photos in clinical environments.

Technical features FALKOLED Illuminator OKKIO HDMI

Model	IL901H
Application	OKKIO HDMI
LED Mount Type	SCREW 1/4-20 UNC
Anti-Polymerization Filter	/

Adjustable Intensity	TOUCHPAD	
Maximum Intensity (LUX)		24000
Color Temperature (K)		5700
CRI	> 80	
Spot Shape	CIRCULAR	
Spot Size at 250 mm	8 CM	
Power Connector	JACK 5.5 x 2.1mm	
Power Supply	12 VDC	
Current Consumption (mA)		500
Cable Length (m)	/	
Battery	/	
Weight (g)		10
Dimensions (mm D x L)	18 x 38	