



Original NSK Torque Wrench Set, 3 Wrenches + 3 Tips NEW

Original NSK torque wrench set including 3 wrenches and 3 G4/G4/G6 tips. Brand new item, availability: 1 piece.

Available

SKU: OFOUT121

Categories: Outlet

Description

Original NSK torque wrench set consisting of 3 torque wrenches and 3 tips G4/G4/G6 for professional dental applications.

New product available in 1 piece only.

Features

- Original NSK set
- 3 torque wrenches included
- 3 tips included: G4/G4/G6
- New product
- Availability: 1 piece only