



Glass cabinet 1000x380x480(h)mm, inox, frosted glasses, 1 glass shelf

Available on request

SKU: DE0097X

Categories: Display Cabinets

Description

The wall-mounted stainless-steel medical cabinet 1000×380×480 H with smoked glass is designed for clinics, dental surgeries and hospital environments that demand hygiene, durability and functional design.

Developed to complement the **Chirunox line**, it provides a **practical and elegant wall-mounted solution** for safe, organised storage.

Constructed entirely from medical-grade stainless steel, it ensures resistance, longevity and easy cleaning, making it ideal for healthcare and clinical use.

The **sliding smoked glass doors** protect stored items while maintaining a modern, professional appearance.

Inside, a **glass shelf** allows efficient organisation of instruments, and the **integrated blue LED lighting** offers soft, even illumination that enhances the perception of cleanliness and calm.

Compact and wall-mounted, with overall dimensions of 1000×380×480 H mm, this stainless-steel medical cabinet is the ideal furnishing solution for clinics and sterile support areas that require hygiene, safety and contemporary design.