



## Glass cabinet 1000x380x480(h)mm, inox, frosted glasses, 1 glass shelf

**Available on request**

**SKU:** DE0097X

**Categories:** Display Cabinets

### Description

---

The **wall-mounted stainless-steel medical cabinet 1000x380x480 H with smoked glass** is designed for **clinics, dental surgeries and hospital environments** that demand **hygiene, durability and functional design**.

Developed to complement the **Chirunox line**, it provides a **practical and elegant wall-mounted solution** for safe, organised storage.

Constructed entirely from **medical-grade stainless steel**, it ensures **resistance, longevity and easy cleaning**, making it ideal for healthcare and clinical use.

The **sliding smoked glass doors** protect stored items while maintaining a modern, professional appearance.

Inside, a **glass shelf** allows efficient organisation of instruments, and the **integrated blue LED lighting** offers soft, even illumination that enhances the perception of cleanliness and calm.

Compact and **wall-mounted**, with overall dimensions of **1000x380x480 H mm**, this **stainless-steel medical cabinet** is the ideal furnishing solution for **clinics and sterile support areas** that require hygiene, safety and contemporary design.