



Glass cabinet 1000x380x480(h)mm,inox,frosted glasses,shelf,dispencer

Available on request

SKU: DE0098X

Categories: Display Cabinets

Description

The stainless-steel medical wall cabinet for disposable products 1000×380×480 H is designed for clinics, dental surgeries and laboratories that demand hygiene, organisation and reliability. Made entirely of medical-grade stainless steel, it ensures durability, corrosion resistance and easy cleaning, ideal for clinical environments.

Specifically created for storing **disposable products** such as gloves, paper or masks, it keeps consumables protected from dust and contamination while maintaining them always within reach.

The **sliding smoked-glass doors** offer a refined appearance and protect the interior from external

The **sliding smoked-glass doors** offer a refined appearance and protect the interior from externa exposure, combining aesthetics and practicality.

Inside, a **central glass shelf** allows efficient organisation of materials, while the **integrated blue LED lighting** provides soft, uniform illumination that enhances the perception of cleanliness and calm. Its **wall-mounted design** and compact dimensions (**1000×380×480 H mm**) help optimise available space, making it ideal for **clinical**, **laboratory and medical areas**.

Combining functionality, hygiene and a clean aesthetic, this **stainless-steel medical wall cabinet for disposable products** is an essential component in any professional healthcare workspace.