



## Glass cabinet 1000x380x480(h)mm,inox,frosted glasses,shelf,dispencer

Wall-mounted stainless steel cabinet for disposable medical supplies 1000x380x480 H mm, with sliding smoked glass doors, internal glass shelf and integrated blue LED light.

**Available on request**

**SKU:** DE0098X

**Categories:** Display Cabinets

### Description

---

The **wall-mounted stainless steel cabinet for disposable medical supplies 1000x380x480 H mm** is a medical furniture unit designed to keep gloves, paper, masks and disposable materials organised, protected and always available in clinics, dental practices and medical laboratories.

Made from **medical-grade stainless steel**, it provides resistance, durability and easy cleaning, essential features in professional environments where hygiene and practical use are a priority.

The wall-mounted structure helps optimise workspace, keeping operating surfaces clear and improving room organisation. The **sliding smoked glass doors** protect the contents from dust and contamination, while giving the cabinet a clean, discreet and professional appearance.

The interior includes a **glass shelf**, useful for separating materials and allowing quick access during daily procedures. This configuration helps organise disposable products efficiently, keeping them ready for use when needed.

The **integrated blue LED light** enhances the cabinet design and helps create a clean, modern and professional environment, conveying a sense of order and safety to the patient.

With dimensions of **1000x380x480 H mm**, the **wall-mounted stainless steel cabinet for disposable medical supplies** is a functional solution for medical environments requiring hygiene, ergonomics and efficient organisation.