


Hanging cabinet 1000x300x500 mm with sliding smoked glass doors



Hanging cabinet with sliding smoked glass doors and two glass shelf, dimensions: 1000 x 300 x 500 mm.

Available on request

SKU: DE0096

Categories: Display Cabinets

Description

The **wall-mounted medical cabinet 1000x300x500 mm with smoked sliding glass** is designed for clinics, dental practices, laboratories and all healthcare environments that require a compact, organised and hygienic storage solution. Its wall-mounted configuration optimises vertical space and keeps work surfaces clear, helping maintain a tidy and efficient clinical setting.

The cabinet features a robust, easy-to-clean structure suitable for modern medical interiors where hygiene and order are essential.

The **smoked sliding glass panels** protect the stored items from dust and external contaminants, while offering a discreet and professional appearance. The smoked finish also helps maintain a visually clean and harmonious environment.

Inside, the cabinet includes **two glass shelves**, allowing instruments and materials to be organised efficiently and accessed quickly during daily procedures. The internal layout ensures clarity, visibility and immediate availability of essential items.

With its clean lines and compact proportions, the wall-mounted medical cabinet with sliding glass integrates seamlessly into a wide range of healthcare facilities, contributing to a modern, orderly and professional atmosphere.