



Hanging cabinet 1000x300x500 mm with sliding smoked glass doors

Hanging cabinet with sliding smoked glass doors and two glass shelf, dimensions: 1000 x 300 x 500 mm.

Available on request

SKU: DE0096

Categories: Sterilization

Description

Hanging cabinet with glass doors.

Tecnomed Italia glass cabinets are indispensable for complete the furniture.

These cabinets offer practicality and functionality: facilitating an efficient and orderly organization.

Products are characterized by a solid and sturdy structure, equipped with sliding frosted glass panel and two glass shelf.

Dimensions: 1000 x 300 x 500 mm.