




Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/hanging-cabinet-1000x300x500-mm-with-sliding-smoked-glass-doors/>



Hanging cabinet 1000x300x500 mm with sliding smoked glass doors

Hanging cabinet with sliding smoked glass doors and two glass shelf, dimensions: 1000 x 300 x 500 mm.

906,51 €

Available on request

SKU: DE0096

Categories: Sterilization

Description

Hanging cabinet with glass doors.

Tecnomed Italia glass cabinets are indispensable for complete the furniture.

These cabinets offer practicality and functionality: facilitating an efficient and orderly organization.

Products are characterized by a solid and sturdy structure, equipped with sliding frosted glass panel and two glass shelf.

Dimensions: 1000 x 300 x 500 mm.