



Hanging cabinet 1000x300x500 mm with sliding smoked glass doors

Hanging cabinet with sliding smoked glass doors and two glass shelf, dimensions: $1000 \times 300 \times 500$ mm.

Available on request

SKU: DE0096

Categories: Display Cabinets, Display Cabinets

Description

Hanging cabinet with glass doors.

Tecnomed Italia glass cabinets are indispensable for compete the furniture.

These cabinets offer praticality and functionality: facilitating an efficent and orderly organization.

Products are characterized by a solid and sturdy structure, equipped with sliding frosted glass panel and two glass shelf.

Dimensions: 1000 x 300 x 500 mm.