

## Helix Device Class 2 test, pack of 100 strips



Helix Device Class 2 test, pack of 100 strips, compatible with all autoclaves

**Available**

**SKU:** DEA4.02

**Categories:** Sterilization Indicators

### Description

---

Helix Device Class 2 test is a physical test for class B autoclaves. It is ideal for assessing the penetration strength of the sterilization agent (steam, hydrogen peroxide, plasma ethylene oxide, formaldehyde) inside hollow bodies.

This test should be carried out in dental offices as there is presence of instruments with cavities of type "A" narrow bodies such as dental handpieces.

Air and other non-condenser gases must be removed from the sterilization chamber before sterilization agent come into contact with the hollow devices to be treated to avoid negatively affecting the outcome of the cycle.

The test is carried out in an empty chamber after the vacuum test and before the normal use of the autoclave.

The Helix Device Class 2 test is compatible with all types of autoclaves and is produced in accordance with current regulations.

- It allows you to immediately and accurately detect the presence of any anomalies during the Helix-test cycle
- The resulting response remains unchanged over time
- Inside a toning indicator, at first glance, indicates whether strip has been used or not.
- 100 strips pack.