



Printed on 30/05/2025

<https://shop.tecnomeditalia.com/en/product/ink-pad-for-labelling-machine-pack-of-5/>

Ink pad for labelling machine, pack of 5


Ink pad for labelling machine, suitable for 2 print lines, diameter 12mm and width 22mm pack of 5

26,01 €

Available

SKU: DEA4.15

Categories: Consumable



Description

Ink pad for labelling machine.

Ink pads with special ink, tightly closed, fitting for labelers with two printing lines.

Expiration date: 3 years

Pack containing 5 pieces