



INOX Hanging cabinet 1500x380x480 mm light colored glass, 1 shelf

Available on request

SKU: DE00975X

Categories: Display Cabinets

Description

The **wall-mounted stainless steel hanging cabinet with clear glass doors 1500x380x480 H mm** is a professional medical furniture unit designed to organise and display materials, instruments and products in clinics, dental practices and medical laboratories.

Made from **medical-grade stainless steel**, it provides resistance, durability and easy cleaning, essential qualities in professional environments where hygiene, order and practicality are important.

The wall-mounted hanging structure helps optimise the workspace, keeping operating surfaces clear and improving room organisation. **Sliding clear glass doors** protect the contents while allowing full visibility of the items stored inside, giving the cabinet a clean, functional and professional appearance.

The interior includes a **glass shelf**, useful for separating materials and making access easier during daily use. This configuration supports efficient organisation and allows the contents to be identified quickly.

With dimensions of **1500x380x480 H mm**, the **wall-mounted stainless steel hanging cabinet with clear glass doors** is a practical and elegant solution for medical environments requiring hygiene, visibility and efficient organisation.