



Inox Luzzani MINIBRIGHT 6F syringe hose,1350mm,for Anthos,SterW,Castel..

Syringe hose, 1350mm, Anthos - Stern Weber - Castellini compatible, 3 functions, complete with inox Minibright lux Luzzani syringe.

Shipping expected in approximately 15 days

SKU: PC714EL

Categories: Syringe

Description

6-function syringe hose compatible with Anthos, Stern Weber, Castellini, length 1350mm, supplied with Luzzani Minibright lux stainless steel syringe for reliable performance.