



## Inox Luzzani MINILIGHT 3F syringe hose,1160mm,for Anthos,SterW,Castel

Syringe hose, 1160mm, Anthos - Stern Weber - Castellini compatible, 3 functions, complete with inox Minilight lux Luzzani syringe

**Shipping expected in approximately 15 days**

**SKU:** PC714C

**Categories:** Syringe

### Description

---

Three-function syringe hose, 1160mm length, compatible with Anthos, Stern Weber, and Castellini units, complete with Minilight stainless steel syringe by Luzzani. High-quality spare part for professional use.