



Syringe hose, 1200mm,3f for KaVo with Minilight Syringe 1057-1058 compat

Syringe hose, 1200mm, 3 functions, KaVo 1057 - 1058 - 1060 - 1062 - 1063 - 1065 - 1066 - 1078 - 1080 compatible, complete with inox Minilight lux Luzzani syringe.

Shipping expected in approximately 15 days

SKU: PC179KN

Categories: Syringe

Description

Syringe hose with 3-function Minilight syringe, 1200mm, KaVo compatible 1057, 1058, 1060, 1062, 1063, 1065, 1066, 1078, 1080 models. Reliable and precise spare part, ideal for fast and secure replacement in dental units.