



Inox Luzzani MINILIGHT 3F syringe hose,1350mm,for Anthos,SterW,Castel

Syringe hose, 1350mm, Anthos - Stern Weber - Castellini compatible, 3 functions, complete with inox Minilight lux Luzzani syringe.

Shipping expected in approximately 15 days

SKU: PC714L

Categories: Syringe

Description

1350mm 3-function syringe hose compatible with Anthos, Stern Weber, Castellini, supplied with Luzzani Minilight stainless steel syringe for high reliability and durability.