



Syringe hose,1500mm,for KaVo 3f,with Minilight inox e70s/e70t/e80t comp.

Syringe hose, 1500mm, 3 functions, KaVo E70S - E70T - E80T compatible, 750ohm, complete with inox Minilight lux Luzzani.

Shipping expected in approximately 15 days

SKU: PC179KN2L

Categories: Syringe

Description

3-function syringe hose compatible with KaVo E70S E70T E80T, 1500mm length, supplied with Minilight Luzzani stainless steel syringe, 750ohm resistance.