


Retrograde root canal preparation endo tip compatible with DTE ED10



Endo tip for retrograde root canal preparation, indicated for finishing and polishing the root canal. Suitable for DTE, Satelec and NSK ultrasonic scalers.

Available

SKU: PG123

Categories: Satelec - DTE - NSK

Description

Endo tip for retrograde root canal preparation. It is indicated for finishing and polishing the root canal during retrograde preparation procedures.

Specifications

- Indicated for canal finishing and polishing
- Compatible with DTE ED10
- Suitable for DTE, Satelec and NSK ultrasonic scalers
- For professional dental use