



Apical surgery endo tip compatible with DTE ED10D, Satelec AS6D,NSK E32D

Diamond-coated endo tip for apical surgery, indicated for retrograde root canal preparation and effective canal finishing. Suitable for DTE, Satelec and NSK ultrasonic scalers.

Available

SKU: PG124

Categories: Satelec - DTE - NSK

Description

Diamond-coated endo tip for apical surgery. It is indicated for retrograde root canal preparation and effective canal finishing.

Specifications

- Diamond-coated endo tip
- Indicated for apical surgery and retrograde preparation
- Compatible with DTE ED10D, Satelec AS6D and NSK E32D
- Suitable for DTE, Satelec and NSK ultrasonic scalers
- For professional dental use