



Endo tip for retrograde root canal preparation, compatible with DTE ED11

Endo tip for retrograde root canal preparation, indicated for controlled finishing and polishing of the root canal. Suitable for DTE, Satelec and NSK ultrasonic scalers.

Available

SKU: PG125

Categories: Satelec - DTE - NSK

Description

Endo tip for scaler, for retrograde root canal preparation, compatible with DTE ED11. It is indicated for controlled finishing and polishing of the root canal.

Specifications

- Compatible with DTE ED11
- Suitable for DTE, Satelec and NSK ultrasonic scalers
- For professional dental use