



Endo preparation tip compatible with DTE ED11D, Satelec AS3D, NSK E31D

E11D Endo tip for scaler, useful for removing very hard materials, SATELEC, DTE, NSK compatible, autoclavable 134°C

Available

SKU: PG126

Categories: Satelec - DTE - NSK

Description

Diamond-coated endo tip for scaler, for retrograde preparation, DTE ED11D, Satelec AS3D and NSK E31D compatible. It is indicated for effective finishing of the root canal during retrograde preparation procedures.

Specifications

- Diamond-coated endo tip
- Compatible with DTE ED11D, Satelec AS3D and NSK E31D
- Suitable for DTE, Satelec and NSK ultrasonic scalers
- For professional dental use