


Broken instrument removal endo tip compatible with DTE ED14, NSK E8



Long endo tip for root canal retreatment, indicated for removing broken instruments inside the root canal. Suitable for DTE, Satelec and NSK ultrasonic scalers.

Available

SKU: PG127

Categories: Satelec - DTE - NSK

Description

Long endo tip for scaler, for root canal retreatment, DTE ED14 and NSK E8 compatible. It is indicated for removing fractured instruments inside the root canal, supporting controlled access along the canal.

Specifications

- Compatible with DTE ED14 and NSK E8
- Suitable for DTE, Satelec and NSK ultrasonic scalers
- For professional dental use